



Piezosurgical tip for scaler, comp. Mectron OT2, Woodpecker US3

Piezosurgical tip for separating the alveolar rim. Compatible Mectron OT2, Woodpecker US3

Available

SKU: PG076M

Categories: Mectron - Woodpecker

Description

Piezosurgery tip for scaler compatible with Mectron OT2 and Woodpecker US3 is designed for alveolar ridge separation in dental surgical procedures.

Specifications

- Compatible with Mectron OT2, Woodpecker US3
- Suitable for use with compatible Mectron and Woodpecker piezoelectric scalers
- For professional dental use