



## Piezosurgical tip for scaler, comp. Mectron OT2, Woodpecker US3

Piezosurgical tip for separating the alveolar rim. Compatible Mectron OT2, Woodpecker US3

**Estimated shipping in 40 days**

**SKU:** PG076M

**Categories:** Mectron - Woodpecker

### Description

---

Piezosurgical tip for scaler, compatible Mectron OT2, Woodpecker US3.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Piezosurgical tip for separating the alveolar rim.

Compatibility Piezosurgical tip for scaler: Mectron OT2, Woodpecker US3.