



Printed on 06/06/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ot5-woodpecker-ul3/>



## Piezosurgical tip for scaler, comp. Mectron OT5, Woodpecker UL3

Diamond disc piezosurgical drill to flatten the vestibular flap bone edge. Compatible Mectron OT5, Woodpecker UL3

**58,41 €**

**Available on request**

**SKU:** PG080M

**Categories:** Piezosurgery inserts

### Description

---

Piezosurgical tip for scaler PG080M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Diamond disc piezosurgical drill to flatten the vestibular flap bone edge.

Compatibility Piezosurgical tip for scaler PG080M.

- Mectron OT5;
- Woodpecker UL3.