



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ot6-woodpecker-us1/>



Piezosurgical tip for scaler, comp. Mectron OT6, Woodpecker US1

Piezosurgical tip for large-section osteotomy as maxillo-facial surgery. Compatible Mectron OT6, Woodpecker US1

70,11 €

Available

SKU: PG074M

Categories: Piezosurgery inserts

Description

Piezosurgical tip for scaler PG074M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Piezosurgical tip for large-section osteotomy as maxillo-facial surgery.

Compatibility Piezosurgical tip for scaler PG074M.

- Mectron OT6;
- Woodpecker US1.