



Piezosurgical tip for scaler, comp. Mectron OT6, Woodpecker US1

Piezosurgical tip for large-section osteotomy as maxillo-facial surgery.
Compatible Mectron OT6, Woodpecker US1

Available

SKU: PG074M

Categories: Mectron - Woodpecker

Description

Piezosurgery tip for scaler compatible with Mectron OT6 and Woodpecker US1, designed for wide-section osteotomy in maxillofacial surgical procedures.

Specifications

- Compatible with Mectron OT6 and Woodpecker US1
- Suitable for use with compatible Mectron and Woodpecker piezoelectric scalers
- For professional dental use