



Printed on 06/06/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ot7-woodpecker-us2/>



Piezosurgical tip for scaler, comp. Mectron OT7, Woodpecker US2

Piezosurgical tip for maxillary and mandibular osteotomy.
Compatible Mectron OT7, Woodpecker US2

70,11 €

Available

SKU: PG075M

Categories: Piezosurgery inserts

Description

Piezosurgical tip for scaler PG075M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Compatibility Piezosurgical tip for scaler PG075M.

- Mectron OT7;
- Woodpecker US2.