



## Piezosurgical tip for scaler, comp. Mectron OT7, Woodpecker US2

Piezosurgical tip for maxillary and mandibular osteotomy. Compatible Mectron OT7, Woodpecker US2

**Available**

**SKU:** PG075M

**Categories:** Mectron - Woodpecker

### Description

---

Piezosurgery tip for scaler compatible with Mectron OT7 and Woodpecker US2, designed for upper maxillary and mandibular osteotomy in dental and maxillofacial surgical procedures.

### Specifications

- Compatible with Mectron OT7 and Woodpecker US2
- Suitable for use with compatible Mectron and Woodpecker piezoelectric scalers
- For professional dental use