

Piezosurgical tip for scaler, comp. Mectron OT8L, Woodpecker US1L



Left-oriented piezosurgical tip for cutting the left monocortical bone.
Comp. Mectron OT8L, Woodpecker US1L

Available

SKU: PG083M

Categories: Mectron - Woodpecker

Description

Piezosurgical tip for scaler, compatible Mectron OT8L, Woodpecker US1L.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Left-oriented piezosurgical tip for cutting the left monocortical bone.

Compatibility Piezosurgical tip for scaler: Mectron OT8L, Woodpecker US1L.