



Piezosurgical tip for scaler, comp. Mectron OT8R, Woodpecker US1R

Right-oriented piezosurgical drill for cutting right monocortical bone.
Compat. Mectron OT8R, Woodpecker US1R

Available

SKU: PG084M

Categories: Mectron - Woodpecker

Description

Piezosurgical tip for scaler compatible with Mectron OT8R and Woodpecker US1R is designed for cutting right monocortical bone.

Specifications

- Compatible with Mectron OT8R and Woodpecker US1R
- Suitable for use with compatible Mectron and Woodpecker piezoelectric scalers
- For professional dental use