



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/piezosurgical-tip-for-scaler-comp-mectron-ot8r-woodpecker-us1r/>



Piezosurgical tip for scaler, comp. Mectron OT8R, Woodpecker US1R

Right-oriented piezosurgical drill for cutting right monocortical bone. Compat. Mectron OT8R, Woodpecker US1R

71,02 €

Available

SKU: PG084M

Categories: Piezosurgery inserts

Description

Piezosurgical tip for scaler PG084M.

Tecnomed Italia stainless steel piezosurgery inserts are reinforced with a titanium nitride coating to improve reliability and ensure superior accuracy in operations on bone.

Traditional tools like burs and saws can easily provoke permanent damage if they come into contact with nerve fibre.

Tecnomed piezosurgery inserts do not cause damage in the event of occasional or accidental contact with nerve fibre, and because bone is not traumatised by cutting, post-operative recovery is faster.

Right-oriented piezosurgical drill for cutting right monocortical bone.

Compatibility

- Mectron OT8R;
- Woodpecker US1R.