



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/pneumatic-veterinary-dental-cart-airo-with-3-instruments/>



## Pneumatic veterinary dental cart Airo with 3 instruments

Airo veterinary dental pneumatic cart with DCI style syringe and 2 M4 cords

**1.996,45 €**

**Available on request**

**SKU:** DE700EP

**Categories:** Veterinary unit

### Description

Pneumatic veterinary dental cart Airo.

Our mobile veterinary dental units are easy to transport and can therefore be used in a number of different veterinary studios.

Their intelligent design guarantees unobstructed operations. Functions are controlled by simple and intuitive hand and foot controls. These units offer all the functions you need to provide effective veterinary dental treatments.

They are made from top quality materials for unrivalled robustness. The spray/water supply from the tank guarantees total operational autonomy. All you need to be up and running is a compressor connection.

The Airo control panel includes: in the upper part the lever selector for dispensing water spray with the relative adjustment tap; in the lower part the 3 lever selectors for choosing the instruments (for example: pneumatic micromotor with contra-angle, high speed handpiece and pneumatic scaler).

The trolley unit come with a progressive foot control. An increase in foot pressure gives an increase in instrument speed.

Technical features Pneumatic veterinary dental cart Airo.

Power Supply	/
Air supply	5-8 bar
Dimensions (w x h x l)	500 x 900 x 500 mm
Weight	15 Kg
DCI style syringe	•
High speed handpiece /pneumatic micromotor cord	• 2

High speed handpiece / pneumatic micromotor cord	Optional
Disc type progressive foot control	•
1lt spray water tank	•
Instrument label	Abs
Warranty	24 months