



Printed on 05/06/2025

<https://shop.tecnomeditalia.com/en/product/pneumatic-veterinary-dental-cart-dentaway-with-3-instruments/>



Pneumatic veterinary dental cart Dentaway with 3 instruments

Pneumatic veterinary dental cart Dentaway with DCI style syringe and 2 M4 cords

1.541,24 €

Available on request

SKU: UV3PN

Categories: Veterinary unit

Description

Pneumatic veterinary dental cart Dentaway.

Our mobile veterinary dental units are easy to transport and can therefore be used in a number of different veterinary studios.

Their intelligent design guarantees unobstructed operations. Functions are controlled by simple and intuitive hand and foot controls. These units offer all the functions you need to provide effective veterinary dental treatments.

They are made from top quality materials for unrivalled robustness. The spray/water supply from the tank guarantees total operational autonomy. All you need to be up and running is a compressor connection.

The front panel houses the water/ spray regulator and a selector for the left or right hand instrument. In the lower part there is the water/ spray lever selector.

The trolley unit come with a progressive foot control. An increase in foot pressure gives an increase in instrument speed.

Technical features Pneumatic veterinary dental cart Dentaway.

Power Supply	/
Air supply	5-8 bar
Dimensions (w x h x l)	400 x 800 x 350 mm
Weight	10 Kg
DCI style syringe	•
High speed handpiece /pneumatic micromotor cord	• 2

High speed handpiece / pneumatic micromotor cord	/
Disc type progressive foot control	•
1lt spray water tank	•
Instrument label	Metal
Warranty	24 months