



Power Pen LED curing light, 1100 mW/cm2, cordless

Power Pen LED curing light, 1100 mW/cm2, cordless, 1 curing programme and 3 curing times

Available

SKU: PP1100W

Categories: cordless dental curing lights

Description

The Power Pen LED curing light 1100 mW/cm2 is the ideal solution for optimal and uniform curing, ensuring high-quality results for all light-curable dental materials.

Description: Designed with cutting-edge technology, the Power Pen combines lightness and precision, allowing you to work efficiently even during extended sessions. Its ergonomic design and user-friendly operation make it particularly suitable for daily use in the dental office.

Strengths/Features:

- Ease of Use: Exposure times can be selected with a single button, with an audible signal every 10 seconds to help you manage curing times accurately.
- Lightweight and Ergonomic: The lightweight design allows prolonged use without causing hand fatigue.
- Autoclavable Tip: Equipped with a black glass fiber optic tip that is autoclavable at 135°C, ensuring maximum hygiene.
- Long-lasting Battery: Memory-free battery, easily replaceable, with an autonomy of 2000 seconds of exposure.
- Curing Modes: Includes "Ramp" mode and three curing times: 10, 20, and 30 seconds.
- Intuitive Indicator Band: The handpiece features an indicator band that changes color, providing immediate feedback on the status of the curing light.
- Stable Base: Comes with a base for secure and stable placement within the dental practice.

Standard equipment:

- Handpiece
- Handpiece base
- 50 disposable dams
- · Battery charger



Technical features Power Pen LED curing light 1100 mW/cm2.

Version	Cordless
Light Source	4 LED modules
Wavelenght	410-490 nm
Light intensity	up to 1100 mW/cm2 +/-5%
Light intensity Hyper setting	3500 mW/cm2 +/-5%
Handpiece dimensions	25 x 208 mm
Lamp tip dimensions	8 mm
Handpiece weight	113 gr.
Power supply	Li-Ion Battery 3,7V (800mAh)
Warranty	12 months