



## Pull-out shelf for Easy Tray, size 415x300x20 mm

Pull-out tray compatible with Easy Tray medical trolleys, designed for the orderly placement of clinical instruments and devices.

**Available on request**

**SKU:** DE1003X.02

**Categories:** Stainless Steel Trays

### Description

---

The **Easy Tray pull-out tray** 415x300x20mm is designed to improve the organisation of instruments and devices on Easy Tray medical trolleys. It allows quick access to equipment during clinical procedures and helps keep the workstation efficient and well organised.

### Key features

- Pull-out tray compatible with Easy Tray medical trolleys
- Suitable for clinical instruments and devices
- Smooth sliding for quick access
- Easy-to-clean surface
- Dimensions: 415x300x20mm