



Pull-out shelf for Easy Tray, size 465X315X20 mm

Pull-out tray compatible with Easy Tray medical trolleys, designed for the orderly placement of clinical instruments and devices.

Available

SKU: DE1003X.01

Categories: Stainless Steel Trays

Description

The **Easy Tray pull-out tray** 465x315x20mm is designed for the organised placement of clinical instruments and devices.

It is suitable for workstations that require practicality, order and efficient management of equipment during clinical activities.

Key features

- Pull-out tray compatible with Easy Tray medical trolleys
- Surface for clinical instruments and devices
- Smooth sliding movement
- Easy to clean and sanitise
- Dimensions: 465x315x20mm