



Injection Tip, 23G 24mm, for Fi-E gutta percha Sealing System

Injection Tip, Type 23G 24mm, for Fi-E Woodpecker Gutta-percha Sealing System.

Available

SKU: E23G-M

Categories: Gutta-Percha Tips

Description

Injection Tip for Fi-E Woodpecker 23G 24mm, for gutta percha Sealing System.

Injection tip for use with heated gutta percha injection pen.

Provides a direct, controlled flow of Gutta Percha and easy reach to the canal.

For use of Backfilling with plasticized GP after the master cone. Type 23G, 24mm.