



## Injection Tip, 23G 28mm, for Fi-E gutta percha Sealing System

Injection Tip, Type 23G 28mm, for Fi-E Woodpecker Gutta-percha Sealing System.

**Available**

**SKU:** E23G-L

**Categories:** Gutta-Percha Tips

### Description

---

Injection Tip for Fi-E Woodpecker 23G 28mm, for gutta percha Sealing System

Injection tip for use with heated gutta percha injection pen.

Provides a direct, controlled flow of Gutta Percha and easy reach to the canal.

For use of Backfilling with plasticized GP after the master cone. Type 23G, 28mm.