



Injection Tip, 25G 28mm, for Fi-E gutta percha Sealing System

Injection Tip for Fi-E Woodpecker 25G 28mm, for Gutta-percha Sealing System.

Available

SKU: E25G-L

Categories: Gutta-Percha Tips

Description

Injection Tip for Fi-E Woodpecker 25G 28mm, for gutta percha Sealing System.

Injection tip for use with heated gutta percha injection pen.

Provides a direct, controlled flow of Gutta Percha and easy reach to the canal.

For use of Backfilling with plasticized GP after the master cone. Type 25G, 28mm.