



Recessed stainless steel washbasin, 500X650X1020h with ster. and valves

Available on request

SKU: WASHCH-3

Categories: Sinks

Description

The **surgical stainless-steel retractable scrub sink** is engineered for **operating theatres, clinics and sterilisation rooms** that require **maximum hygiene, safety and functionality**.

Made entirely of **stainless steel**, it combines **durability, corrosion resistance and easy cleaning**, ensuring compliance with healthcare hygiene standards.

The **retractable design** allows the sink to be safely concealed when not in use, protecting the surface and optimising workspace hygiene.

Its **foot-operated solenoid valve** provides **hands-free water activation**, improving aseptic workflow and preventing cross-contamination.

This precise, contactless control system is ideal for **pre-surgical handwashing procedures** in sterile areas.

An integrated **UVC sterilisation system** delivers continuous antibacterial and antiviral action on internal surfaces, maintaining high hygiene levels throughout clinical use.

With dimensions of **500x650x1020 mm**, the **deep stainless-steel bowl** features **rounded edges and smooth welded surfaces**, eliminating residue build-up and simplifying sanitation.

Compact, hygienic and robust, this **surgical stainless-steel retractable scrub sink** offers a **professional-grade solution** for hospitals and surgical departments focused on hygiene, safety and efficiency.