



Recessed stainless steel washbasin,800X650X1020h with ster. and valves

Available on request

SKU: WASHCH-380

Categories: Sinks

Description

STERILA is a **surgical stainless-steel retractable scrub sink 800mm**, designed for **operating theatres, clinics and sterilisation rooms** where **maximum hygiene and space efficiency** are essential. Made entirely of **stainless steel**, it ensures **strength, durability and easy cleaning**, meeting the strict hygiene standards required in healthcare facilities.

The **retractable design** allows the sink to be concealed when not in use, protecting it from contaminants and maintaining a clean, uncluttered workspace.

The **foot-operated solenoid valve** enables **hands-free activation**, ensuring precise control of the water flow and reducing the risk of cross-contamination. The **UVC sterilisation system** provides effective antibacterial and antiviral action, keeping the sink hygienically safe between uses.

With overall dimensions of **800x650x1020 H**, the **deep stainless-steel bowl** offers generous space for surgical handwashing.

Its **800 mm width** ensures comfort and accessibility, while the **smooth welded surfaces and rounded corners** prevent residue build-up and simplify sanitation.

Compact, reliable and hygienic, **STERILA**, a **surgical stainless-steel retractable scrub sink 800mm** represents a **professional solution for hospitals and surgical departments** that demand safety, ergonomics and efficiency.