



Recessed stainless steel washbasin,800X650X1020h with ster. and valves

Available on request

SKU: WASHCH-380 Categories: Sinks

Description

STERILA is a surgical stainless-steel retractable scrub sink 800mm, designed for operating theatres, clinics and sterilisation rooms where maximum hygiene and space efficiency are essential. Made entirely of stainless steel, it ensures strength, durability and easy cleaning, meeting the strict hygiene standards required in healthcare facilities.

The **retractable design** allows the sink to be concealed when not in use, protecting it from contaminants and maintaining a clean, uncluttered workspace.

The **foot-operated solenoid valve** enables **hands-free activation**, ensuring precise control of the water flow and reducing the risk of cross-contamination. The **UVC sterilisation system** provides effective antibacterial and antiviral action, keeping the sink hygienically safe between uses.

With overall dimensions of **800×650×1020 H**, the **deep stainless-steel bowl** offers generous space for surgical handwashing.

Its **800 mm width** ensures comfort and accessibility, while the **smooth welded surfaces and rounded corners** prevent residue build-up and simplify sanitation.

Compact, reliable and hygienic, STERILA, a surgical stainless-steel retractable scrub sink 800mm represents a professional solution for hospitals and surgical departments that demand safety, ergonomics and efficiency.