

Rotor, comp., W&H SYNEA TA-96, ceramic bearing, Myonic, torlon



Available

SKU: WT33CC

Categories: Synea-Alegra

Description

Compatible rotor for W&H Synea TA-96 Mini and W&H Alegra TE-96BC high speed handpiece, supplied with Myonic torlon ceramic bearings.

Handpiece compatibility

- W&H Synea TA-96 Mini
- W&H Alegra TE-96BC

Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Marked with serial number
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards