



Rotor compatible Castellini EVO standard, ceramic phenolic bearings

Rotor for dental high speed handpiece Ccompatible Castellini EVO standard, replacement part supplied with phenolic ceramic bearings.

Available on request

SKU: HCL843CC

Categories: Castellini

Description

Compatible Castellini EVO standard rotor is a replacement part for dental high speed handpieces, supplied with phenolic ceramic bearings.

Compatible handpieces

- Castellini EVO standard

Technical features

- Phenolic ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur retention >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349