


## Rotor, compatible Castellini Titanium Gold HP2, Myonic ceramic bearings



Compatible rotor for Castellini Hi Power Ceramic 2, Titanium Gold 2 / L high speed handpieces, for standard head models, supplied with Myonic torlon ceramic bearings.

**Available on request**

**SKU:** HCC862C

**Categories:** Castellini

### Description

---

Compatible rotor for Castellini Hi-Power 2 ceramics – Titanium Gold/2L high speed handpiece, replacement part for standard head models, supplied with Myonic torlon ceramic bearings.

### Compatibility

- Castellini Hi Power Ceramic 2
- Castellini Titanium Gold 2

### Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur retention force >40 Newton
- 12-month warranty
- Complies with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349