


Rotor, compatible Castellini Titanium Gold HP2, Myonic ceramic bearings



Compatible rotor for Castellini Hi Power Ceramic 2, Titanium Gold 2 / L high speed handpieces, for standard head models, supplied with Myonic torlon ceramic bearings.

Available

SKU: HCC862C

Categories: Castellini

Description

Compatible rotor for Castellini Hi-Power 2 ceramics – Titanium Gold/2L high speed handpiece, replacement part for standard head models, supplied with Myonic torlon ceramic bearings.

Compatibility

- Castellini Hi Power Ceramic 2
- Castellini Titanium Gold 2

Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur retention force >40 Newton
- 12-month warranty
- Complies with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349