

Rotor compatible head ENDOMATE NSK /ENDOSMART MM18EN



Complete rotor compatible with NSK Endo-Mate MP-F16R 16:1 and MP-F20R 20:1 heads.

Available

SKU: 521TM015

Categories: NSK

Description

NSK Endo-Mate MP-F16R MP-F20R head rotor complete compatible, designed to restore precise performance in endodontic low speed handpiece heads.

Compatibility

- NSK Endo-Mate MP-F16R 16:1
- NSK Endo-Mate MP-F20R 20:1

Features

- Precisely balanced
- 12-month warranty
- Compliance with EN ISO 13485 standards