

## Rotor compatible head KAVO CL3



Compatible rotor for KaVo SURGmatic S201L C3 CL3 rotor, with steel bearings.

**Estimated shipping in 30 days**

**SKU:** K10N4490

**Categories:** Kavo

### Description

---

Compatible KaVo SURGmatic S201L C3 CL3 rotor, engineered spare part for restoring proper operation of dental instruments.

### Compatibility

- KaVo SURGmatic S201L
- KaVo C3
- KaVo CL3
- Ivory SL631C

### Features

- Steel bearings
- Perfectly balanced
- 12-month warranty
- Compliant with EN ISO 13485 standards