



Rotor, compatible KaVo Gentleforce 7000B/C, Barden steel torlon bearings

Compatible rotor for KaVo Gentleforce 7000B and 7000C dental high speed handpiece, supplied with Myonic torlon steel bearings.

Available on request

SKU: HKT700T

Categories: GF 7000

Description

Complete compatible rotor for KaVo Gentleforce 7000B and 7000C, replacement part for high speed handpiece, supplied with Myonic torlon steel bearing.

Handpiece compatibility

- KaVo Gentleforce 7000B
- KaVo Gentleforce 7000C

Original KaVo rotor reference code: **1.000.4772**

Technical features

- Myonic torlon steel bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349