



## Rotor, compatible Midwest Tradition,without bearings

Compatible rotor for Midwest Tradition L and PB dental high speed handpiece, supplied without bearings.

**Available**

**SKU:** SC101

**Categories:** Midwest

### Description

---

Complete compatible rotor for Midwest Tradition L and PB, spare part for dental high speed handpiece supplied without bearings.

#### Handpiece compatibility

- Midwest Tradition L
- Midwest Tradition PB

#### Technical features

- Dynamically balanced rotor
- Marked with serial number
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349