



## Rotor compatible Sirona T1 Boost K.SN > 800K, ceramic torlon bearings

Compatible rotor for Sirona T1 Boost K high speed handpiece, s.n. > 800K, supplied with ceramic torlon bearings.

**Available on request**

**SKU:** R464SCH

**Categories:** Boost

### Description

---

Compatible rotor for **Sirona T1 Boost K, s.n. > 800K**, replacement part for dental high speed handpieces supplied **with ceramic torlon bearings**.

#### Handpiece compatibility

- Sirona T1 Boost K
- Sirona T1/T2/T3 Boost
- SN > 800.000 / s.n. > 800K

#### Technical features

- ceramic torlon bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Manufactured in compliance with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349