



Rotor compatible Sirona T1 Boost. Sn > 800K. Ceramic torlon bearings

Compatible rotor for Sirona T2/T3 Boost high speed handpiece, supplied with Myonic torlon ceramic bearings, for units with SN over 800000.

Available

SKU: R464SK

Categories: Boost

Description

The **compatible rotor for Sirona T2 / T3 Boost high speed handpiece dental turbine** is a spare part supplied with **Myonic torlon ceramic bearings**.

Handpiece compatibility

- Sirona T1 Boost*
- Sirona T2 Boost*
- Sirona T3 Boost*

*with SN over 800000

Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Identified with serial number
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards