



Rotor compatible Sirona T2 CONTROL, Myonic ceramic bearings

Compatible rotor for Sirona T1 T2 Control dental high speed handpiece, dental turbine spare part supplied with Myonic torlon ceramic bearings.

Available

SKU: R054TK

Categories: Control

Description

Compatible rotor for Sirona T1 / T2 Control dental high speed handpiece, dental turbine spare part supplied with Myonic torlon ceramic bearings and is designed for versions with 3-notch push button.

Handpiece compatibility

- Sirona T1 Control with serial number above 600000 and below 800000
- Sirona T2 Control with serial number above 600000 and below 800000

Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards