


## Rotor, compatible Sirona T2/T3 Racer, Myonic ceramic torlon bearings



Compatible rotor for Sirona T2 / T3 Racer high speed handpiece, Sirona Siroboost S/T, B.A. International BA695L, B.A. International BA755L and Fona 8080L SR, with Myonic torlon ceramic bearings.

**Available on request**

**SKU:** R046RC

**Categories:** Racer

### Description

---

The complete compatible rotor for Sirona T2 / T3 Racer is a replacement part for a dental high speed handpiece, supplied with Myonic torlon ceramic bearings.

#### Handpiece compatibility

- Sirona T2 Racer
- Sirona T3 Racer
- Sirona Siroboost S/T
- B.A. International BA695L
- B.A. International BA755L
- Fona 8080L SR

Original SIRONA spare part number: **5868307**

#### Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur holding force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1 and EN ISO 5349 standards