



Rotor, compatible Sirona TS1-2 4000S, Myonic ceramic torlon bearings

Compatible rotor for Sirona TS1, TS2 and 4000S high speed handpiece, dental turbine spare part supplied with Myonic torlon ceramic bearings.

Available on request

SKU: RC024S

Categories: 4000 S-MS

Description

Compatible rotor for Sirona TS1 TS2 4000S high speed handpiece, spare part for dental turbine supplied with Myonic torlon ceramic bearings.

Handpiece compatibility

- Sirona TS1
- Sirona TS2
- Sirona 4000S

Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur holding force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349