



## Rotor, compatible W&H Synea TA-97, Myonic ceramic torlon bearings

**Available on request**

**SKU:** WTA97LCM

**Categories:** Synea-Alegra

### Description

---

**Compatible rotor for W&H Synea TA-97 LED and Alegra TE-97 high speed handpiece**, spare part for dental turbine supplied with Myonic torlon ceramic bearings.

#### Handpiece compatibility

- W&H Synea TA-97
- W&H Synea TA-97 LED
- W&H Alegra TE-97 old model

#### Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards