

## Rotor for head,original Nsk TI95 (revision), completo



Complete rotor for NSK Ti-Max TI95 Rev.A low speed handpiece, suitable for instruments with serial number above A:G83, with Myonic torlon ceramic bearings.

**Available**

**SKU:** NS311N

**Categories:** NSK

### Description

---

Complete rotor for NSK Ti-Max TI95 Rev.A low speed handpiece, compatible replacement for contra-angles with serial number above A:G83.

### Compatibility

- NSK Ti-Max TI95 Rev.A with serial number above A:G83

### Features

- Myonic torlon ceramic bearings
- Perfectly balanced
- 12-month warranty
- Complies with EN ISO 13485 standards