

Rotor for head,original Nsk TI95 (revision), complete



Complete rotor for NSK Ti-Max TI95 Rev.A low speed handpiece, suitable for instruments with serial number above A:G83, with Myonic torlon ceramic bearings.

Available

SKU: NS311N

Categories: NSK

Description

Complete rotor for NSK Ti-Max TI95 Rev.A low speed handpiece, compatible replacement for contra-angles with serial number above A:G83.

Compatibility

- NSK Ti-Max TI95 Rev.A with serial number above A:G83

Features

- Myonic torlon ceramic bearings
- Perfectly balanced
- 12-month warranty
- Complies with EN ISO 13485 standards