



## Rotor Ivory, complete, for S504 4 spray, from the SN 2501001

Compatible complete rotor for Ivory S504M high speed handpiece, 4 spray, from SN 2501001, supplied with torlon ceramic bearings.

**Estimated shipping in 15 days**

**SKU:** TIV500.3

**Categories:** Ivory

### Description

---

Complete rotor for Ivory S504M 4 spray from SN 2501001, spare part for dental high speed handpiece, supplied with torlon ceramic bearings.

#### Handpiece compatibility

- Ivory S504M 4 spray with serial number ? 2501001

#### Technical features

- Supplied with torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards