


Rotor, MK-dent small 1 spray comp., Myonic ceramic torl. bearing, Tecnomed



Compatible rotor for MK-dent HC4014 HC6012 1 spray dental high speed handpiece, supplied with Myonic torlon ceramic bearings.

Available

SKU: TC212T

Categories: Classic Line

Description

Compatible rotor for MK-dent HC4014 HC6012 1 spray dental high speed handpiece is supplied with Myonic torlon ceramic bearing.

Turbine compatibility

- HC4012K with SN up to E
- HC4012M with SN up to E
- HC4012S with SN up to E
- HC6012K with SN up to E
- HC6012S with SN up to E

MK-dent rotor reference code: **TC2012**.

Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur holding force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349