

Rotor, original TKD Thera L-R-E



Original rotor for TKD Thera L-R-E dental high speed handpiece, supplied with torlon steel bearings.

Shipping expected in approximately 15 days

SKU: 184.90

Categories: TKD

Description

Original rotor for TKD Thera L – R – E, spare part for dental high speed handpiece, supplied with torlon steel bearings.

Handpiece compatibility

- TKD Thera L
- TKD Thera R
- TKD Thera E

Referal TKD rotor code: **184.90**

Technical features

- Supplied with torlon steel bearings
- Dynamically balanced rotor
- Serial number identification
- Bur holding force >40 Newton
- 12-month warranty
- Manufactured in compliance with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349