



Rotor, original TKD Titan, Titan Lux, Thera, Thera Lux

Compatible original rotor for TKD Titan, Titan Lux, Thera and Thera Lux high speed handpiece supplied with torlon steel bearings.

Available

SKU: 180.90

Categories: TKD

Description

Original rotor for TKD Titan, Titan Lux, Thera and Thera Lux high speed handpiece, dental turbine spare part supplied with torlon steel bearings.

Compatible handpieces

- TKD Titan
- TKD Titan Lux
- TKD Thera
- TKD Thera Lux

TKD rotor code: **180.90**

Technical features

- Torlon steel bearings
- Dynamically balanced rotor
- Serial number identification
- Bur holding force >40 Newton
- 12-month warranty
- Manufactured in compliance with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349