



Rotor, compatible MK-dent Prime Line standard, Myonic ceramic bearing

Compatible rotor for **MK-dent Prime Line standard dental high speed handpiece**, supplied with torlon ceramic bearings.

Available on request

SKU: TC241T

Categories: Prime Line

Description

Complete compatible rotor for **MK-dent Prime Line standard**, spare part for dental high speed handpiece, supplied with torlon ceramic bearings.

Handpiece compatibility

- MK-dent HC8021B
- MK-dent HC8021K
- MK-dent HC9021NT

MK-dent rotor reference code: **TC2041**

Technical features

- torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards