


## Rotor, compatible W&H 698-798-898-190 B/C, Myonic ceramic torlon bearings



Compatible rotor for W&H 698, 798, 898, 195, 190B and 190C high speed handpiece, dental turbine spare part supplied with Myonic torlon ceramic bearings.

**Available on request**

**SKU:** WTX98CC

**Categories:** 3xx-6xx-7xx-8xx

### Description

---

**Compatible rotor for W&H 698, 798, 898, 195, 190B and 190C high speed handpiece,** dental turbine spare part supplied with Myonic torlon ceramic bearings.

#### Turbine compatibility

- W&H 698
- W&H 798
- W&H 898
- W&H 195
- W&H 190B
- W&H 190C

#### Technical features

- Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349