



## Rotor, compatible W&H Synea TA-98/Alegra, 2 grips, Myonic ceramic bearings

Compatible rotor for W&H Synea TA-98 and TA-98LC high speed handpiece with 2-notch push button, supplied with Myonic torlon ceramic bearings.

**Available**

**SKU:** WT05CC

**Categories:** Synea-Alegra

### Description

---

**Compatible rotor for W&H Synea TA-98 and TA-98LC high speed handpiece**, replacement part for dental turbine supplied with Myonic torlon ceramic bearings.

#### Turbine compatibility

- W&H Synea TA-98 with 2-notch push button
- W&H Synea TA-98LC with 2-notch push button

#### Technical features

- Supplied with Myonic torlon ceramic bearings
- Dynamically balanced rotor
- Marked with serial number
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards