

## Rotor, compatible W&H Synea TA-98/Alegria, Myonic steel torl. bear., 2 grips



Compatible rotor for W&H Synea TA-98 and TA-98LC high speed handpiece with 2-notch push button, supplied with Myonic torlon steel bearings.

**Available**

**SKU:** WT05S

**Categories:** Synea-Alegria

### Description

---

**Compatible rotor for W&H Synea TA-98 and TA-98LC high speed handpiece**, replacement part for dental turbine supplied with Myonic torlon steel bearings.

#### Turbine compatibility

- W&H Synea TA-98
- W&H Synea TA-98LC

#### Technical features

- Supplied with Myonic torlon steel bearings
- Dynamically balanced rotor
- Marked with serial number
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards