


Rotor, compatible W&H Synea TA-98/Alegra, 2 grips, Myonic steel bearings



Compatible rotor for W&H Synea TA-98 and TA-98LC high speed handpiece with 2-notch push button, supplied with Myonic torlon steel bearings.

Available

SKU: WT05S

Categories: Synea-Alegra

Description

Compatible rotor for W&H Synea TA-98 and TA-98LC high speed handpiece, replacement part for dental turbine supplied with Myonic torlon steel bearings.

Turbine compatibility

- W&H Synea TA-98
- W&H Synea TA-98LC

Technical features

- Supplied with Myonic torlon steel bearings
- Dynamically balanced rotor
- Marked with serial number
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards