



## Rotor, MK-dent Classic Line 4 spraystandard compatible, ceramic bearing

Compatible rotor for MK-dent Classic Line standard 4 spray dental high speed handpiece, supplied with torlon ceramic bearings, compatible with HE17, HE21K, HE21KL and HB21KL.

**Available on request**

**SKU:** TC251T

**Categories:** Classic Line

### Description

---

Compatible rotor for MK-dent Classic Line standard 4 spray, replacement rotor for standard head dental high speed handpiece, supplied with torlon ceramic bearings.

#### Handpiece compatibility

- MK-dent Classic Line HE17
- MK-dent Classic Line HE21K
- MK-dent Classic Line HE21KL
- MK-dent Classic Line HB21KL

MK-dent rotor reference code: **TC1034**

#### Technical features

- Torlon ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur holding force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349