



Rotor,SIRONA T3-T4 RACER SN>600.000 compatible,ceramic phenolic bearing

Compatible rotor for Sirona T3 T4 Racer high speed handpiece with Myonic phenolic ceramic bearings.

Available on request

SKU: R461SK

Categories: Racer

Description

Compatible rotor for Sirona T3 / T4 Racer high speed handpiece is a replacement part for dental turbine supplied with Myonic phenolic ceramic bearings.

Turbine compatibility

- Sirona T3 Racer with serial number between 600000 and 700000
- Sirona T4 Racer with serial number between 600000 and 700000

Technical features

- Myonic phenolic ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards