



## Rotor, SIRONA T3-T4 RACER SN > 600.000 compatible, ceramic phenolic bearing

Compatible rotor for Sirona T3 T4 Racer high speed handpiece with Myonic phenolic ceramic bearings.

**Available**

**SKU:** R461SK

**Categories:** Racer

### Description

---

**Compatible rotor for Sirona T3 / T4 Racer high speed handpiece** is a replacement part for dental turbine supplied with Myonic phenolic ceramic bearings.

#### Turbine compatibility

- Sirona T3 Racer with serial number between 600000 and 700000
- Sirona T4 Racer with serial number between 600000 and 700000

#### Technical features

- Myonic phenolic ceramic bearings
- Dynamically balanced rotor
- Serial number identification
- Bur gripping force >40 Newton
- 12-month warranty
- Compliant with EN ISO 13485, EN ISO 10993-1, EN ISO 17665-1, EN ISO 5349 standards