



Printed on 31/05/2025

<https://shop.tecnomeditalia.com/en/product/round-seat-various-colours-chrome-plated-base-tecnomed/>



## ROUND seat, various colours, chrome-plated base, Tecnomed

ROUND seat, various colours, chrome-plated base, Tecnomed height adjustment

**592,42 €**

**Available**

**SKU:** DE8.08X

**Categories:** Handy Line

### Description

---

Round seat.

An appropriate stool is an important investment for the comfort and health of the dental professional.

Beauty professional dentist chair is designed for complete comfort even after many hours of work.

The ergonomic structure helps you to maintain correct posture while distributing weight in a uniform manner.

The antibacterial and water-repellent imitation leather upholstery makes these chairs extremely easy to disinfect.

The range offers a choice of thirty different colours for a perfect match with the finish and materials of your dental chair.

Top quality imitation leather upholstery ensures a soft, anatomical seating position with uniform weight distribution and no painful compressions.

Round seat feature adjustable pneumatic height adjustment and backrest depth height and angle adjustment.